

Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA

Prepared for:

Nick Jeffery
Aquatic Science Biologist
Fisheries and Oceans Canada
Nick.Jeffery@dfo-mpo.gc.ca

RPT037-OP01

July 15, 2022



Centre for Environmental Genomics Applications



**Centre for Environmental
Genomics Applications**


eDNAtec Inc.

103-14 International Pl.
St. John's NL A1A 0R6
CANADA

T: +1 (709) 576-2550

services@ednatec.com
www.cegacanada.com

This page is intentionally left blank

	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

The correct citation for this report is:

Center for Environmental Genomics Applications, Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA, CEGA RPT037-OP01, 2022-07-15.

Project Team:


Mehrdad Hajibabaei, MH (Director)

Nicole Fahner, NF (Associate Director)

Eman Elbakry, EE (Research Technologist)

Muneesh Kaushal, MK (Bioinformatician)

Jasmine Yeung, JY (Junior Ecologist)

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

REVISION HISTORY

VERSION	NAME	COMPANY	DATE OF ISSUE	COMMENTS
R1	Lesley Berghuis	CEGA	2022-07-15	First release

DISTRIBUTION LIST

NAME	COMPANY	NUMBER OF COPIES / FORMAT
Mehrdad Hajibabaei	CEGA	1 / electronic [pdf]
Nicole Fahner	CEGA	1 / electronic [pdf]
Nick Jeffery	DFO	1/electronic [pdf]



 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15


TABLE OF CONTENTS

1	Executive Summary.....	1
2	Background	2
3	Methods.....	3
3.1	Sample Submission	3
3.2	Sample Submission	3
3.3	Data Analysis.....	3
4	Results.....	6
5	References	7

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15


LIST OF TABLES

Table 1. List of water eDNA filters received for analysis and the DNA extract concentration.	4
Table 2. List of water filter field blanks received for analysis and the DNA extract concentrations. BQL indicates that the DNA in that sample was below detection limit.	5
Table 3. List of DNA extracts received for analysis and their concentrations. BQL indicates that the DNA in that sample was below detection limit.	5

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

1 Executive Summary

Fisheries and Oceans Canada (DFO) is using environmental genomics for harbour monitoring in Halifax, NS. In 2021, DFO collected environmental water samples from harbours in Halifax and submitted them to CEGA for laboratory analysis. 100 Sterivex filters (96 environmental samples and 4 field blanks) and six DNA extracts were submitted and received for analysis. DNA was extracted from all filters and two DNA markers (12S MiFishU and 16S from McInnes et al. 2017) were amplified and then sequenced on an Illumina NovaSeq 6000 instrument. DNA extraction and sequencing was successful for all samples. Raw sequencing results were demultiplexed and trimmed and provided to DFO for further bioinformatics analysis.

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15


2 Background

Environmental genomics is a novel approach to biodiversity characterization that does not require collection of whole biological specimens but instead relies on recovery and analysis of DNA from the physical environment in which they live (e.g., water, soil, sediment, etc.). This environmental DNA (eDNA) is released from organisms through various mechanisms including cell shedding and the excretion of various bodily fluids and feces¹⁻³. DNA is isolated from samples, biodiversity informative regions of DNA are amplified through polymerase chain reaction (PCR), and then these regions are sequenced using high-throughput genomic sequencing platforms.

There are several advantages to environmental genomics compared to the conventional approach to biodiversity characterization. It is non-invasive and does not require the capture of a whole biological specimen. Very small samples of water (~250 mL) or sediment (~5 g) are required for a biodiversity assessment because, unlike conventional surveys, organisms at all trophic levels from algae to large mammals can be detected from the same environmental sample. The technique is highly sensitive, which means rare or elusive species, including endangered and invasive species, can be detected through their DNA in addition to the more common species. It is cost effective with time savings during field sampling and sample analysis. Furthermore, high-throughput DNA sequencing technology allows for simultaneous analysis of a large set of samples.

Because environmental genomics is a novel approach to biodiversity characterization, it is not widely applied yet. However, in recent years, there has been a steady increase in its application to various environmental characterization efforts⁴⁻⁶. While the sensitivity of environmental genomics is one of the many appeals of this technology, it is also one of its potential limitations as contaminants of biological nature can also be detected within the environmental sample. To manage this, many precautions are taken during the acquisition, preservation, and processing of samples to limit the introduction of contaminants (e.g. collection of field blanks as a negative control).

For this project, Fisheries and Oceans Canada (DFO) is using environmental genomics for harbour monitoring, particularly of fish biodiversity, in Halifax, Canada.

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

3 Methods

3.1 Sample Submission

1L environmental water samples were collected and filtered by DFO from harbours in Halifax in August 2021 and submitted to CEGA for lab analysis. CEGA received 100 Sterivex filters (96 environmental samples and 4 field blanks) preserved in Longmire's buffer and six DNA extracts. All samples were received frozen and were stored at -80°C until processing.

3.2 Sample Processing

DNA was extracted from filters for each sample and field blank and quantified for all samples, including those received as DNA extracts (Tables 1-3). Two DNA markers for fish biodiversity, 12S MiFishU and the 16S marker from McInnes et al. (2017)⁷, were amplified in all samples and blanks. Additional negative controls were added during extraction and PCR stages, which were carried through subsequent stages. All samples, including field and lab negative controls, were sequenced on an Illumina NovaSeq 6000 instrument using a 2 x 150 cycle SP kit and a minimum target depth of 1,000,000 reads per sample per marker.

3.3 Data Analysis

Minimal processing was done to produce the fastq files from the raw sequencer data. Demultiplexing was performed using Illumina's BCL Convert software (version 3.10.5). Primers were then trimmed from the forward and reverse reads using cutadapt v4.1. The maximum error rate was set to 0.07 with a minimum overlap equal to 13 bases. Read pairs in which the primer was missing from either the forward or reverse read were discarded.


 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

Table 1. List of water eDNA filters received for analysis and the DNA extract concentration.

Sample Name	[DNA] (pg/ μ L)	Sample Name	[DNA] (pg/ μ L)	Sample Name	[DNA] (pg/ μ L)
485525	621	485618	2984	485661	454
485526	1650	485619	1277	485662	1517
485532	383	485620	1609	485663	598
485533	1995	485621	1251	485664	756
485539	741	485622	2455	485666	388
485540	1213	485624	721	485667	1896
485543	707	485625	793	485668	1006
485544	1049	485626	1254	485669	1801
485545	807	485627	1745	485670	331
485546	1893	485628	454	485671	2353
485547	523	485629	1375	485673	919
485548	2082	485631	781	485674	2762
485563	1792	485632	3257	485675	948
485564	1882	485633	905	485676	2548
485568	334	485634	3719	485677	603
485569	2161	485635	845	485678	4198
485573	822	485636	4361	485680	1285
485574	1821	485638	758	485681	2504
485589	1135	485639	5165	485682	635
485590	1644	485640	741	485683	2820
485594	1928	485641	4376	485684	1511
485595	1063	485642	509	485685	2750
485599	3961	485643	3666		
485600	1714	485645	1098		
485603	500	485646	2665		
485604	317	485647	1133		
485605	572	485648	2817		
485606	1722	485649	1375		
485607	325	485650	3687		
485608	1922	485652	1176		
485610	454	485653	3852		
485611	1456	485654	1358		
485612	902	485655	4887		
485613	1716	485656	1699		
485614	781	485657	4691		
485615	1705	485659	615		
485617	1867	485660	911		



 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

Table 2. List of water filter field blanks received for analysis and the DNA extract concentrations. BQL indicates that the DNA in that sample was below detection limit.

Sample Name	[DNA] (pg/ μ L)
485609	22
485658	BQL
485686	BQL
485549	BQL

Table 3. List of DNA extracts received for analysis and their concentrations. BQL indicates that the DNA in that sample was below detection limit.

Sample Name	[DNA] (pg/ μ L)
DFOGUL21_41102	51
DFOGUL21_41103	91
DFOGUL21_41104	265
DFOGUL21_41105	BQL
DFOGUL21_41106	1138
DFOGUL21_41107	279

 Centre for Environmental Genomics Applications	Monitoring Harbour Fish Communities in Halifax, Nova Scotia using Environmental DNA	
	RPT037-OP01	2022-07-15

4 Results

The mean DNA concentration recovered from the submitted filters was 1664 pg/μL (range: 317 – 5165 pg/μL, Table 1) and all submitted field blanks had DNA concentrations below quantifiable levels, except for one which had a DNA concentration of 22 pg/μL (Table 2). The mean DNA concentration of the submitted DNA extracts was 365 pg/μL (range: 51 – 1138 pg/μL with one sample below the quantification limit (Table 3)). There was an average of 1,321,008 sequencing reads (range: 39,969 – 5,860,704 reads) per sample per marker, excluding lab and field negative controls.

In summary, DNA was extracted and sequenced successfully for all samples. Files with raw sequence data for further bioinformatics analysis will be provided to DFO.

5 References

1. Baird, D. J. & Hajibabaei, M. Biomonitoring 2.0: a new paradigm in ecosystem assessment made possible by next-generation DNA sequencing. *Mol. Ecol.* **21**, 2039–2044 (2012).
2. Taberlet, P., Coissac, E., Hajibabaei, M. & Rieseberg, L. H. Environmental DNA. *Mol. Ecol.* **21**, 1789–1793 (2012).
3. Thomsen, P. F. *et al.* Detection of a diverse marine fish fauna using environmental DNA from seawater samples. *PLoS ONE* **7**, e41732 (2012).
4. Hajibabaei, M., Shokralla, S., Zhou, X., Singer, G. A. C. & Baird, D. J. Environmental barcoding: a next-generation sequencing approach for biomonitoring applications using river benthos. *PLoS ONE* **6**, e17497 (2011).
5. Leray, M. & Knowlton, N. Censusing marine eukaryotic diversity in the twenty-first century. *Philos. Trans. R. Soc. B Biol. Sci.* **371**, 20150331 (2016).
6. Gibson, J. F. *et al.* Large-Scale Biomonitoring of Remote and Threatened Ecosystems via High-Throughput Sequencing. *PLOS ONE* **10**, e0138432 (2015).
7. McInnes, J. C. *et al.* DNA Metabarcoding as a Marine Conservation and Management Tool: A Circumpolar Examination of Fishery Discards in the Diet of Threatened Albatrosses. *Front. Mar. Sci.* **4**, 277 (2017).